

Кронштейны

Технические характеристики

По вопросам продаж и поддержки обращайтесь:

Алматы (727)345-47-04
Ангарск (3955)60-70-56
Архангельск (8182)63-90-72
Астрахань (8512)99-46-04
Барнаул (3852)73-04-60
Белгород (4722)40-23-64
Благовещенск (4162)22-76-07
Брянск (4832)59-03-52
Владивосток (423)249-28-31
Владикавказ (8672)28-90-48
Владимир (4922)49-43-18
Волгоград (844)278-03-48
Вологда (8172)26-41-59
Воронеж (473)204-51-73
Екатеринбург (343)384-55-89

Иваново (4932)77-34-06
Ижевск (3412)26-03-58
Иркутск (395)279-98-46
Казань (843)206-01-48
Калининград (4012)72-03-81
Калуга (4842)92-23-67
Кемерово (3842)65-04-62
Киров (8332)68-02-04
Коломна (4966)23-41-49
Кострома (4942)77-07-48
Краснодар (861)203-40-90
Красноярск (391)204-63-61
Курск (4712)77-13-04
Курган (3522)50-90-47
Липецк (4742)52-20-81

Магнитогорск (3519)55-03-13
Москва (495)268-04-70
Мурманск (8152)59-64-93
Набережные Челны (8552)20-53-41
Нижний Новгород (831)429-08-12
Новокузнецк (3843)20-46-81
Ноябрьск (3496)41-32-12
Новосибирск (383)227-86-73
Омск (3812)21-46-40
Орел (4862)44-53-42
Оренбург (3532)37-68-04
Пенза (8412)22-31-16
Петрозаводск (8142)55-98-37
Псков (8112)59-10-37
Пермь (342)205-81-47

Ростов-на-Дону (863)308-18-15
Рязань (4912)46-61-64
Самара (846)206-03-16
Санкт-Петербург (812)309-46-40
Саратов (845)249-38-78
Севастополь (8692)22-31-93
Саранск (8342)22-96-24
Симферополь (3652)67-13-56
Смоленск (4812)29-41-54
Сочи (862)225-72-31
Ставрополь (8652)20-65-13
Сургут (3462)77-98-35
Сыктывкар (8212)25-95-17
Тамбов (4752)50-40-97
Тверь (4822)63-31-35

Тольятти (8482)63-91-07
Томск (3822)98-41-53
Тула (4872)33-79-87
Тюмень (3452)66-21-18
Ульяновск (8422)24-23-59
Улан-Удэ (3012)59-97-51
Уфа (347)229-48-12
Хабаровск (4212)92-98-04
Чебоксары (8352)28-53-07
Челябинск (351)202-03-61
Череповец (8202)49-02-64
Чита (3022)38-34-83
Якутск (4112)23-90-97
Ярославль (4852)69-52-93

Россия +7(495)268-04-70

Казахстан +7(727)345-47-04

Беларусь +(375)257-127-884

Узбекистан +998(71)205-18-59

Киргизия +996(312)96-26-47

эл.почта: mqc@nt-rt.ru || сайт: <https://milliq.nt-rt.ru/>



WBPRPKL01

Pack Holder for PrePak L Pretreatment Pack

Securely supports the PrePak L Pretreatment pack

packagingpkg of 1 unit
H × W × D795 mm (31.3 in.) × 234 mm (9.2 in.) × 192 mm (7.6 in.)
compatibilityfor use with Elix®
.....for use with Milli-Q® IX 7003
.....for use with Milli-Q® IX 7005
.....for use with Milli-Q® IX 7010
.....for use with Milli-Q® IX 7015
.....for use with RiOs™
shipped inambient

Description

General description

The PrePak L wall mounting bracket is a water purification accessory designed to support the PrePak L pretreatment pack for a leak-free and secure fitting.

The unit includes wall bracket, fittings, tubing, and pressure gauges required for the installation and operation of the PrePak L pretreatment pack that protects Elix® and RiOs™ water purification systems.

Features and Benefits

Tight and leak-free fitting of PrePak L Pretreatment pack

Legal Information

ELIX is a registered trademark of Merck KGaA, Darmstadt, Germany

Milli-Q is a registered trademark of Merck KGaA, Darmstadt, Germany

RIOS is a trademark of Merck KGaA, Darmstadt, Germany



ZFPACKSP2 Pack Support 2 Tube

with 2 gauges, For secure, leak-free fit of the pretreatment pack

packagingpkg of 1 unit
compatibility.....for use with AFS®
shipped inambient

Description

General description

The Pack Support (2 tube) is a water purification accessory that comes with two pressure gauges. It provides a leak-free fit of replacement cartridge for MILLI-RO™ Compact Water Systems. The accessory is needed for installation and operation of PrePak 1 pretreatment pack.

Legal Information

AFS is a registered trademark of Merck KGaA, Darmstadt, Germany
MILLI-RO is a trademark of Merck KGaA, Darmstadt, Germany



ZFPACKSP0 Pack Support 2 Tube

without gauges, For secure, leak-free fit of the pretreatment pack

packagingpkg of 1 unit
compatibility.....for use with Direct-Q® UV
.....for use with Direct-Q®

Description

General description


The Pack Support 2 tube is a water purification accessory that provides a leak-free fit of pretreatment pack.

Application

Wall mounting bracket for installation and operation of pretreatment pack

Legal Information

DIRECT-Q is a registered trademark of Merck KGaA, Darmstadt, Germany



TANKFIX01

Wall Mounting Bracket for 30 / 60 L Polyethylene Storage Tanks

Save lab bench space

feature Save lab bench space
packaging pkg of 1 unit
compatibility for use with Direct-Q® UV
..... for use with Elix® Essential UV
..... for use with Elix® Essential Water Purification Systems
..... for use with Milli-Q® Direct Water Purification Systems
..... for use with RiOs™ Essential Water Purification Systems
..... for use with RiOs™ Water Purification System
shipped in ambient

Description

General description

Save lab bench space
Water purification accessory

Application

Save lab bench space

Legal Information

DIRECT-Q is a registered trademark of Merck KGaA, Darmstadt, Germany
ELIX is a registered trademark of Merck KGaA, Darmstadt, Germany
Milli-Q is a registered trademark of Merck KGaA, Darmstadt, Germany
RIOS is a trademark of Merck KGaA, Darmstadt, Germany



WMBWASH1

Wall Mounting Bracket

Washer Distribution Kit installation accessory for use with Elix® Advantage and RiOs™, Saves lab bench space

packagingpkg of 1 unit
compatibility.....for use with Elix® Advantage
.....for use with RiOs™

Description

General description

This accessory is used for wall mounting the washer distribution kit on either the right or left side of the polyethylene (PE) reservoir.

Legal Information

ELIX is a registered trademark of Merck KGaA, Darmstadt, Germany
RIOS is a trademark of Merck KGaA, Darmstadt, Germany



MPKMNTKIT

Wall Mounting Bracket

Millipak® Wall Mounting Bracket saves lab bench space

packagingpkg of 1 unit
manufacturer/tradenameMillipak®
W × H × D.....3.5 in. (89 mm) × 3.1 in. (79 mm) × 1.25 in. (32 mm)
compatibility.....for use with AFS®
shipped inambient

Description

General description

The wall mounting bracket saves lab space by easily installing the Millipak® filter on the wall.

Application

Water purification system accessory for use with AFS® D water purification system.

Legal Information

AFS is a registered trademark of Merck KGaA, Darmstadt, Germany
MILLIPAK is a registered trademark of Merck KGaA, Darmstadt, Germany



MBPKMNKIT
Wall Mounting Bracket

Wall mounting bracket Kit for installing BioPak® C ultrafiltration cartridge

sterility	0.2 µm filtered
feature	removal of bacterial by-products, such as alkaline
.....	phosphatase (ALP) and endotoxins from pure water, 1.5 L/min at 1 bar and 20°C
flow rate	
packaging	pkg of 1 unit
manufacturer/tradename	BioPak® C
parameter	5 °C temp. range (35 °C)
resistivity.....	>10.0 MΩ-cm (at 25 degree celcius)
impurities.....	<10 cfu/mL bacterial count
total organic carbon (TOC) residue	<500 ppb
pore size.....	0.22 µm pore size
fitting.....	1/4 in. inlet connection (Gaz (Straight) Luer Slip F)
.....	8 mm outlet connection (Teflen)
compatibility.....	for use with AFS®
shipped in	dry ice

Description

General description

Wall mounting bracket kit is a water purification accessory with all the necessary components that provide convenient installation of the BioPak® C filter unit.
microbiological filter, does not remove ALP (as expected. Ultrafiltration is efficient for ALP removal)

Application


For use with AFS®, AFS® E and AFS® D Water Purification Systems.

Components

- Metal bracket with adhesive strips
- ¼ in. GazF - 8 mm stem connector
- 8 mm–6 mm T reducer
- 8 mm–6 mm elbow (Qty: 2)
- Installation instruction insert

Legal Information

AFS is a registered trademark of Merck KGaA, Darmstadt, Germany
BIOPAK is a registered trademark of Merck KGaA, Darmstadt, Germany



WMBWASHA1

Wall Mounting Bracket

Installation accessory for use with Milli-Q® IQ and Milli-Q® IX systems

feature	Save lab bench space
packaging	pkg of 1 unit
compatibility.....	for use with Milli-Q® IX 7003
.....	for use with Milli-Q® IX 7005
.....	for use with Milli-Q® IX 7010
.....	for use with Milli-Q® IX 7015
.....	for use with Milli-Q® IQ 7003
.....	for use with Milli-Q® IQ 7005
.....	for use with Milli-Q® IQ 7010
.....	for use with Milli-Q® IQ 7015

Description

General description

The mounting bracket easily installs a Washer Distribution Kit high up on a wall, leaving benchtops free of clutter for optimized lab space and a more efficient working environment.

This bracket is compatible with:

Washer Distribution Kit 230 V (right)

Washer Distribution Kit 230 V (left)

Washer Distribution Kit 115 V (right)

Washer Distribution Kit 115 V (left)

Features and Benefits

Easy to install.

Saves lab bench space.

Legal Information

Milli-Q is a registered trademark of Merck KGaA, Darmstadt, Germany



WMBSMT002
Wall Mounting Bracket

For installation of Milli-Q Advantage/Integral/Reference/Direct, Simplicity, Synergy and Elix Reference systems

feature	Saves lab bench space
packaging	pkg of 1 unit
L × W × H	352 mm × 291 mm × 154 mm
compatibility.....	for use with Elix® Advantage
.....	for use with Elix® Reference
.....	for use with Milli-Q® Advantage A10
.....	for use with Milli-Q® Direct
.....	for use with Milli-Q® Integral
.....	for use with Milli-Q® Reference
.....	for use with Simplicity® UV
.....	for use with Simplicity®
.....	for use with Synergy® UV
.....	for use with Synergy®
shipped in	ambient

Description

General description

Attaching the water purification system to the wall is now feasible with our wall mounting bracket. This easy-to-install water purification accessory for Milli-Q® systems guarantees to save precious lab bench space. The wall mounting bracket supports Milli-Q® Advantage/Integral/Reference/Direct, Simplicity®, Synergy® and Elix® Reference water purification systems.

Features and Benefits

Save lab bench space

Legal Information

ELIX is a registered trademark of Merck KGaA, Darmstadt, Germany
Milli-Q is a registered trademark of Merck KGaA, Darmstadt, Germany
SIMPLICITY is a registered trademark of Merck KGaA, Darmstadt, Germany
SYNERGY is a registered trademark of Merck KGaA, Darmstadt, Germany



SYSTFIX01
Wall Mounting Bracket

Installation accessory for RiOs™ and Elix® water purification systems

packaging	pkg of 1 unit
compatibility.....	for use with AFS®
.....	for use with Elix® Essential UV
.....	for use with Elix® Essential
.....	for use with RiOs™ Essential

Description

General description

The mounting bracket easily installs a water purification unit high up on a wall, leaving benchtops free of clutter for optimized lab space and a more efficient working environment.
This bracket is compatible with RiOs™ Essential and Elix® Essential water purification Systems.

Application

Save lab bench space.

Features and Benefits

Easy to install accessory.
Saves lab bench space.

Legal Information

AFS is a registered trademark of Merck KGaA, Darmstadt, Germany
ELIX is a registered trademark of Merck KGaA, Darmstadt, Germany
RIOS is a trademark of Merck KGaA, Darmstadt, Germany



M0PMNTKIT Wall Mounting Bracket

Saves lab bench space

AC/DC input	100 V / 240 V AC, 50-60 Hz
sterility	0.22 µm filtered
packaging	pkg of 1 unit
manufacturer/tradename	Opticap®
parameter	>5 L/min flow rate at 2 bar
.....	.28 kg operating weight
.....	.5 °C temp. range (35 °C)
resistivity.....	>10.0 MΩ-cm (at 25 degree celcius)
H × W × D.....	585 mm × 268 mm × 426 mm
impurities.....	<10 cfu/mL Bacteria level
total organic carbon (TOC) residue	<500 ppb
pore size.....	0.22 µm pore size
conductivity.....	100-2000 µS/cm at 25 °C
compatibility.....	for use with AFS®
shipped in	ambient

Description

General description

Wall mounting bracket enables easy installation of an Opticap® onto the wall, freeing your bench for crucial research work.

Application

For use with AFS® and AFS® E Water Purification Systems.

Legal Information

AFS is a registered trademark of Merck KGaA, Darmstadt, Germany
OPTICAP is a registered trademark of Merck KGaA, Darmstadt, Germany



WMBSMT001
Wall Mounting Bracket

For installation of RiOs™ 3, RiOs-DI® 3 and Direct-Q® 3 systems

feature	Saves lab bench space
packaging	pkg of 1 unit
compatibility.....	for use with Direct-Q® UV
.....	for use with Direct-Q®
.....	for use with RiOs-DI® UV
.....	for use with RiOs-DI®
.....	for use with RiOs-DI® Clinical
.....	for use with RiOs™
shipped in	ambient

Description

General description

Wall mounting bracket helps you to easily install the water purification system, anywhere you need pure water, on the benchtop or high on the wall, freeing your bench for crucial research work.

Application

Wall Mounting Bracket with Tap supports RiOs™, RiOs-DI® and Direct-Q® Systems.

Features and Benefits

- Increase access to water from your bench.
- Saves more space.

Legal Information

DIRECT-Q is a registered trademark of Merck KGaA, Darmstadt, Germany
RIOS is a trademark of Merck KGaA, Darmstadt, Germany
RIOS-DI is a registered trademark of Merck KGaA, Darmstadt, Germany



WMBQP0D01
Wall Mounting Bracket

Installation accessory for use with POD dispensers

packaging	pkg of 1 unit
H × W × D	31 cm × 15 cm × 44.1 cm
compatibility	for use with E-POD®
.....	for use with Elix® Advantage
.....	for use with Milli-Q® IQ Element
.....	for use with Milli-Q® IX 7003
.....	for use with Milli-Q® IX 7005
.....	for use with Milli-Q® IX 7010
.....	for use with Milli-Q® IX 7015
.....	for use with Q-POD® Element

Description

General description

This wall mounting accessory is designed for easy and uncluttered installation of remote water dispensing units. The bracket includes a rotating arm for further flexibility.
For use with Q-POD® and E-POD® dispensers, and the Q-POD® Element unit.

Features and Benefits

Saves lab bench space.
Easy to install.

Legal Information

E-POD is a registered trademark of Merck KGaA, Darmstadt, Germany
ELIX is a registered trademark of Merck KGaA, Darmstadt, Germany
Milli-Q is a registered trademark of Merck KGaA, Darmstadt, Germany
Q-POD is a registered trademark of Merck KGaA, Darmstadt, Germany



TANKFIXA1 Wall Mounting Bracket

Installation accessory for 25 or 50 L water storage unit, For use with Milli-Q® IQ and Milli-Q® IX systems

packaging	pkg of 1 unit
compatibility	for use with Milli-Q® IX 7003
.....	for use with Milli-Q® IX 7005
.....	for use with Milli-Q® IX 7010
.....	for use with Milli-Q® IX 7015
.....	for use with Milli-Q® IQ 7003
.....	for use with Milli-Q® IQ 7005
.....	for use with Milli-Q® IQ 7010
.....	for use with Milli-Q® IQ 7015

Description

General description

The mounting bracket easily installs a 25 or 50 L water storage unit high up on a wall, leaving benchtops free of clutter for optimized lab space and a more efficient working environment

This bracket installs Milli-Q® IQ 7003/05/10/15 and Milli-Q® IX 7003/05/10/15 storage tanks.

Application

For flexible installation and saving lab bench space.

Features and Benefits

Easy to install accessory.

Saves lab bench space.

Legal Information

ELIX is a registered trademark of Merck KGaA, Darmstadt, Germany

Milli-Q is a registered trademark of Merck KGaA, Darmstadt, Germany



SYSTFIXA1
Wall Mounting Bracket

Installation accessory for use with Milli-Q® IQ, Milli-Q® IX and Milli-Q® EQ water purification systems

packaging	pkg of 1 unit
compatibility.....	for use with Milli-Q® EQ 7008
.....	for use with Milli-Q® EQ 7016
.....	for use with Milli-Q® IX 7003
.....	for use with Milli-Q® IX 7005
.....	for use with Milli-Q® IX 7010
.....	for use with Milli-Q® IX 7015
.....	for use with Milli-Q® IQ 7000
.....	for use with Milli-Q® IQ 7003
.....	for use with Milli-Q® IQ 7005
.....	for use with Milli-Q® IQ 7010
.....	for use with Milli-Q® IQ 7015

Description

General description

The mounting bracket easily installs a water purification unit high up on a wall, leaving benchtops free of clutter for optimized lab space and a more efficient working environment.


This bracket is compatible with Milli-Q® IQ 7000, Milli-Q® IQ 7003/05/10/15, Milli-Q® IX 7003/05/10/15, and Milli-Q® EQ 7000/08/16 water purification systems.

Features and Benefits

- Easy to install accessory.
- Saves lab bench space.

Legal Information

Milli-Q is a registered trademark of Merck KGaA, Darmstadt, Germany



WMEQ0DKT

Milli-Q® EQ Wall Mounting Kit

Wall mounting kit for the HMI touchscreen of Milli-Q® EQ 7000/08/16 water purification systems

compatibility.....for use with Milli-Q® EQ 7000
.....for use with Milli-Q® EQ 7008
.....for use with Milli-Q® EQ 7016

Description

General description

The wall mounting kit is a water system accessory that includes all necessary components for convenient installation of the HMI touchscreen of Milli-Q® EQ 7 series water purification systems. The touchscreen display can be wall-mounted up to 3 m from the system, saving lab space.

This kit for the HMI optimizes the installation flexibility of the system. If your lab changes, the system set-up can adapt to your lab.

Features and Benefits

Ease of installation.

Saves lab bench space.

Legal Information

Milli-Q is a registered trademark of Merck KGaA, Darmstadt, Germany



WMEQ0RKT

Milli-Q® EQ Wall Mounting Kit

Wall mounting kit for the Q-POD® of Milli-Q® EQ series water purification systems.

packagingpkg of 1 unit
compatibility.....for use with Milli-Q® EQ 7000
.....for use with Milli-Q® EQ 7008
.....for use with Milli-Q® EQ 7016

Description

General description

The wall mounting kit is a water system accessory that includes all necessary components for convenient installation of the Q-POD® dispenser of Milli-Q® EQ 7000/08/16 water purification systems. The dispenser can be wall-mounted up to 3 m from the system, saving space.

This kit optimizes the installation flexibility of the system. If your lab changes, the system set-up can adapt to your lab.

Features and Benefits

Saves lab space.
Ease of installation.

Legal Information

Milli-Q is a registered trademark of Merck KGaA, Darmstadt, Germany
Q-POD is a registered trademark of Merck KGaA, Darmstadt, Germany



По вопросам продаж и поддержки обращайтесь:				
Алматы (727)345-47-04 Ангарск (3955)60-70-56 Архангельск (8182)63-90-72 Астрахань (8512)99-46-04 Барнаул (3852)73-04-60 Белгород (4722)40-23-64 Благовещенск (4162)22-76-07 Брянск (4832)59-03-52 Владивосток (423)249-28-31 Владикавказ (8672)28-90-48 Владимир (4922)49-43-18 Волгоград (844)278-03-48 Вологда (8172)26-41-59 Воронеж (473)204-51-73 Екатеринбург (343)384-55-89	Иваново (4932)77-34-06 Ижевск (3412)26-03-58 Иркутск (395)279-98-46 Казань (843)206-01-48 Калининград (4012)72-03-81 Калуга (4842)92-23-67 Кемерово (3842)65-04-62 Киров (8332)68-02-04 Коломна (4966)23-41-49 Кострома (4942)77-07-48 Краснодар (861)203-40-90 Красноярск (391)204-63-61 Курск (4712)77-13-04 Курган (3522)50-90-47 Липецк (4742)52-20-81	Магнитогорск (3519)55-03-13 Москва (495)268-04-70 Мурманск (8152)59-64-93 Набережные Челны (8552)20-53-41 Нижний Новгород (831)429-08-12 Новокузнецк (3843)20-46-81 Ноябрьск (3496)41-32-12 Новосибирск (383)227-86-73 Омск (3812)21-46-40 Орел (4862)44-53-42 Оренбург (3532)37-68-04 Пенза (8412)22-31-16 Петрозаводск (8142)55-98-37 Псков (8112)59-10-37 Пермь (342)205-81-47	Ростов-на-Дону (863)308-18-15 Рязань (4912)46-61-64 Самара (846)206-03-16 Санкт-Петербург (812)309-46-40 Саратов (845)249-38-78 Севастополь (8692)22-31-93 Саранск (8342)22-96-24 Симферополь (3652)67-13-56 Смоленск (4812)29-41-54 Сочи (862)225-72-31 Ставрополь (8652)20-65-13 Сургут (3462)77-98-35 Сыктывкар (8212)25-95-17 Тамбов (4752)50-40-97 Тверь (4822)63-31-35	Тольятти (8482)63-91-07 Томск (3822)98-41-53 Тула (4872)33-79-87 Тюмень (3452)66-21-18 Ульяновск (8422)24-23-59 Улан-Удэ (3012)59-97-51 Уфа (347)229-48-12 Хабаровск (4212)92-98-04 Чебоксары (8352)28-53-07 Челябинск (351)202-03-61 Череповец (8202)49-02-64 Чита (3022)38-34-83 Якутск (4112)23-90-97 Ярославль (4852)69-52-93
Россия +7(495)268-04-70	Казахстан +7(727)345-47-04	Беларусь +(375)257-127-884	Узбекистан +998(71)205-18-59	Киргизия +996(312)96-26-47

эл.почта: mqc@nt-rt.ru || сайт: <https://milliq.nt-rt.ru/>